

# Monte Carlo

## 1- Demonstration of Casino

Open the Casino program

Chose File->New (Edit Layer window)

Click Add Layer

Double click on Undefined (right side of the red circle)

In the composition edit box add C (capitalize)

Click Done button (Layer Chemical Composition window)

Click Next button (Edit Layers window)

Put 5 in the energy field

Put 2000 in the number of electrons to simulate

Click Next button (Microscope and Simulation Proprieties window)

Click Next button (Distributions window)

Click Next button (Options window)

Click Finish button (Select Physics Model window)

Click Yes button (New Simulation Request window)

Wait for the end of the simulation

Click on automatic rescale (View->Auto-Adjust Scale menu)

BSE Coefficient

$$\eta = \underline{\hspace{2cm}}$$

You can note the max depth of the electron and radius

$$D_{\max} = \underline{\hspace{2cm}}$$

$$R_{\max} = \underline{\hspace{2cm}}$$

Show Energy by Position distribution (tree view at the left of the screen)

$$D_{\max} = \underline{\hspace{2cm}}$$

$$R_{\max} = \underline{\hspace{2cm}}$$

Show the phirhoz curves and note the depth and  $\phi(0)$ .

$$D_{\max} = \underline{\hspace{2cm}}$$

$$\phi(0) = \underline{\hspace{2cm}}$$

Now change the specimen for Ag and the energy to 4 keV

Menu Simulation->Modify Sample

Double click on C (right side of the red circle)

In the composition field add Ag

In the name field add Ag

Click Done button

Menu Simulation->Set Up Microscope

Change 5 keV to 4 keV

Click Ok button

Menu Simulation->Run or the green “play” button

BSE Coefficient

$$\eta = \underline{\hspace{2cm}}$$

You can note the max depth of the electron and radius

$$D_{\max} = \underline{\hspace{2cm}}$$

$$R_{\max} = \underline{\hspace{2cm}}$$

Show Energy by Position distribution (tree view at the left of the screen)

$$D_{\max} = \underline{\hspace{2cm}}$$

$$R_{\max} = \underline{\hspace{2cm}}$$

Show the phirhoz curves and note the depth and  $\phi(0)$ .

$$D_{\max} = \underline{\hspace{2cm}}$$

$$\phi(0) = \underline{\hspace{2cm}}$$

## 2- Incident Energy Variation

Simulate for Ag specimen with energies between 5 to 30 keV with a 5 keV step

Menu Simulation->Set Up Microscope

Click Simulate Multiple Energies

Change 4 keV to 5 keV in Start, 30 keV in End and Step to 5 keV

Click Ok button

Menu Simulation->Run or the green “play” button

<b><i>Energy (keV)</i></b>	<b><i><math>\eta</math></i></b>	<b><i><math>D_{\max} e^-</math></i></b>	<b><i><math>R_{\max} e^-</math></i></b>	<b><i><math>D_{\max} x\text{-ray}</math></i></b>	<b><i><math>\phi(0)</math></i></b>
		nm	nm	nm	
5					
10					
15					
20					
25					
30					

Simulate for C specimen with energies between 5 to 30 keV with a 5 keV step

Menu Simulation->Modify Sample

Double click on Ag (right side of the red circle)

In the composition field add C

In the name field add C

Click Done button

Menu Simulation->Run or the green “play” button

<i>Energy (keV)</i>	$\eta$	<i>D<sub>max</sub> e<sup>-</sup></i>	<i>R<sub>max</sub> e<sup>-</sup></i>	<i>D<sub>max</sub> x-ray</i>	$\phi(\theta)$
		nm	nm	nm	
5					
10					
15					
20					
25					
30					

### 3- Atomic Number Variation

Simulate for C, Al, Cu, Ag, and Au specimen with energy of 10 keV at 0°

<i>Element</i>	$\eta$	<i>D<sub>max</sub> e<sup>-</sup></i>	<i>R<sub>max</sub> e<sup>-</sup></i>	<i>D<sub>max</sub> x-ray</i>	$\phi(\theta)$
C					
Al					
Cu					
Ag					
Au					

### 4- Incident Electron Angle

Simulate for Al specimen with energies between 5 keV at 0°, 20°, and 60°

Menu Simulation->Modify Sample

Double click on Au (right side of the red circle)

In the composition field add Al

In the name field add A1

Click Done button

Menu Simulation->Set Up Microscope

Change the Tilt of the specimen to 20°

Click Ok button

Menu Simulation->Run or the green “play” button

<b><i>Tilt (°)</i></b>	<b><i><math>\eta</math></i></b>	<b><i><math>D_{max} e^-</math></i></b>	<b><i><math>R_{max} e^-</math></i></b>	<b><i><math>D_{max} x\text{-ray}</math></i></b>	<b><i><math>\phi(\theta)</math></i></b>
0					
20					
60					

Look also at the Energy by Position distribution

# Microanalysis

## **1- Demonstration of Win X-Ray**

Open the Win X-Ray program

Chose File->New

Change Incident Energy Start to 10 keV

Set the number of electrons to 500

Check X-Ray Compute and Compute Characteristic, and Compute Background

Click Next button (Option Simulation window)

Click on Set Element for All Regions button (Option Specimen window)

Change the atomic number Z to 29, or Symbol to Cu

Click Ok button (Option element for the region 1 window)

Click Next button (Option Specimen window)

Click Next button (Option X-Ray window)

Click Next button (Advanced Option window)

Click Next button (Physics Model window)

Click Next button (Result General Option window)

Click Next button (Result Trajectory Option window)

Click Finish button (Result Distribution Option window)

Wait for the end of the simulation

Click on + to expand the tree view

Click on the trajectory +

Click on the interaction volume and adjust the new windows as you wish

Move the mouse over the picture to find the max depth and radial

Change the X-Z Plane to X-Y Plane

Click on the X-Ray +

Click on the Intensity to shows the intensity generated (I Gen) and emitted (I Emi) (with specimen absorption) for each line in number of photons.

Click on Spectrum to show the complete spectrum with the effect of the detector

Click on the Y Axis Log

Check and uncheck Background, Characteristic, and Total

Move the mouse to read energy and intensity at the cursor.

Without Y log scale, zoom on the K-Lines (8 keV) by left click + move down and right with the mouse to select the zoom region.

To reset the view left-click + move up and left with the mouse.

Now click on PRZ Curves


Click on Ka1 of the Cu, you have to click on + Cu first.

Click on La of the Cu, you should see the effect of the specimen absorption

Click on lines, Cu

Click on Ka1 and Ka2

Click on La and Lb1



## ***2- L lines from Ag***

Select in the menu Option->Option Specimen


Click on Set Element for All Regions button

Change the atomic number Z to 47, or Symbol to Ag

Click Ok button

Click on Simulation->Run or the green “play” button

Look at the Spectrum for the L-lines



### 3- M lines from W

Select in the menu Option->Option Specimen


Click on Set Element for All Regions button

Change the atomic number Z to 74, or Symbol to W

Click Ok button

Click on Simulation->Run or the green “play” button

Look at the Spectrum for the W-lines



#### 4- k-ratio Calculation

Open WinX-Ray

Menu Option->Option Simulation

Set the energy to 15 keV

Check Compute X-Ray

Set the number of electrons to 500

Simulate for pure Au, note the intensity of  $L_{\alpha}$  line

Simulate for pure Cu, note the intensity of  $K_{\alpha 1}$  line

Menu Option->Option Specimen

Click Set Element for All Regions

Set the number of elements at 2 (click on up arrow)

Set element 1 to Au

Set element 2 to Cu

Set Weight Fraction to 0.8 and 0.2 respectively

Click Ok button

Click on Simulation->Run or the green "play" button

Note the intensity for Au  $L_{\alpha}$  and Cu  $K_{\alpha 1}$

Change the Weight Fraction to 0.6 and 0.4 and re simulate

<i>Concentration of Au (%)</i>	<i><math>I L_{\alpha} Au</math></i>	<i>k ratio Au</i>	<i><math>I K_{\alpha 1} Cu</math></i>	<i>k ratio Cu</i>
100				
80				
60				
40				
20				
0				

If k ratio is equal to 0.62 gold and 0.35 for Cu,

what it is the concentration of this alloy? \_\_\_\_\_

#### 5- Thin Film on Substrate

Open Casino

Open Simulation->Set Up Microscope

Put 20 keV of Energy and 2000 for the Number of electrons to simulate



Open Simulation->Modify Sample

Select Multi-Layer

Click Add Layer button

Click to change it to Cu

Click Add Layer button

Click to change it to Si

Change the thickness of Cu to 10 nm

Click Done

Click on Simulation->Run or the green “play” button

Pure Cu Intensity = \_\_\_\_\_

Pure Si Intensity = \_\_\_\_\_

<i>Thickness of Cu (nm)</i>	<i>I K Cu</i>	<i>k ratio Cu</i>	<i>I K Si</i>	<i>k ratio Si</i>
10				
50				
100				
500				
1000				
1500				

## **6- Particles in Matrix**

Open Casino

Open Simulation->Set Up Microscope

Click Simulate Multiple Energies

Put to 1 keV in Start, 10 keV in End and Step to 1 keV

Put 2000 for the Number of electrons to simulate

Open Simulation->Modify Sample

Select Multi-Layer

Click Add Layer button

Click to change it to Ti

Click Add Layer button

Click to change it to Si

Click Add Layer button

Click to change it to Ti

Change the thickness of Ti to 100 nm (in the first layer)

Change the thickness of Si to 100 nm

Click Done

Click on Simulation->Run or the green “play” button

Simulate pure Ti

Simulate pure Si

<b><i>Incident Energy (keV)</i></b>	<b><i>I L Ti</i></b>	<b><i>Pure Ti</i></b>	<b><i>k ratio Ti</i></b>
1			
2			
3			
4			
5			
6			
7			
8			
9			
10			

<b><i>Incident Energy (keV)</i></b>	<b><i>I K Si</i></b>	<b><i>Pure Si</i></b>	<b><i>k ratio Si</i></b>
1			
2			
3			
4			
5			
6			
7			
8			
9			
10			

# How to Solve Problems with Simulations

## **1- Vertical Interface C Au**

Open Casino

Open Simulation->Modify Sample

Select Grain Boundary

Click Add Layer button

Click to change it to C

Click Add Layer button

Click to change it to Au

Click Done button

Open Simulation->Set Up Microscope

Put 20 keV for Energy and 2000 for the Number of electrons to simulate

Click on Simulation->Run or the green “play” button

## **2- Scanning a Vertical Interface**

Open Simulation->Set Up Microscope

Check Scan the beam to create an image

From: -400 nm, To: 400 nm, step: 200 nm

Click on Simulation->Run or the green “play” button

<b>Position</b>	$\eta$
-400	
-300	
-200	
-100	
0	
100	
200	
300	
400	

Contrast = \_\_\_\_\_

## **3- Scanning a Vertical Interface part 2**

Open Casino

Open Simulation->Modify Sample

Select Grain Boundary

Click Add Layer button

Click to change it to GaAs

Click Add Layer button

Click to change it to AlAs

Click Add Layer button

Click to change it to GaAs

Change the thickness of AlAs to 5 nm

Click Done

Open Simulation->Set Up Microscope

Put 15 keV for Energy and 2000 for the Number of electrons to simulate

Put a beam diameter of 2 nm

Check Scan the beam to create an image

From: -10 nm, To: 10 nm, step: 2 nm

Click on Simulation->Run or the green "play" button

<b><i>Position</i></b>	$\eta$
-10	
-8	
-6	
-4	
-2	
0	
2	
4	
6	
8	
10	

Contrast = \_\_\_\_\_

Change the energy to 1 keV with a beam diameter of 5 nm and re simulate

<i>Position</i>	$\eta$
-10	
-8	
-6	
-4	
-2	
0	
2	
4	
6	
8	
10	

Contrast = \_\_\_\_\_